

ORDERING INFORMATION

Cat. №	Product
HC-T2-F	Thrombodynamics Analyser System T2-F (fibrin registration only)
HC-T2-T	Thrombodynamics Analyser System T2-T (fibrin and thrombin registration)
HC-TDX-10	Reagents kit for 10 fibrin measurements
HC-PLS-10	Reagents kit for 10 fibrin and thrombin measurements

The product was launched in the European market as a RUO device at the beginning of 2014.

Thrombodynamics Analyser System is sold in Russia from the end of 2012. Thrombodynamics assay has already been approved for clinical use in Russian Federation.



We help clinicians to succeed, providing clear and reliable estimation of patients coagulation state, thus giving an ability to prevent thrombotic or bleeding complications.

Partners

Switzerland
ENDOTELL AG
info@endotell.ch
www.endotell.ch

Benelux
NODIA
info@nodia.com
www.nodia.com

UK and Channel Islands
Quadrantech Diagnostics Ltd
quadrantech@btinternet.com
www.quadrantech.co.uk

Austria and Germany
CoaChrom Diagnostica GmbH
info@coachrom.com
www.coachrom.com

Don't find your local distributor?
Do not hesitate to contact us:
mail@hemacore.com

www.thrombodynamics.com
www.hemacore.com/en/




HEMOSTASIS & THROMBOSIS BEYOND BIOCHEMISTRY NEW GLOBAL VIEW



THROMBODYNAMICS ASSAY

SPATIAL DYNAMICS OF
FIBRIN CLOT GROWTH



THROMBODYNAMICS-4D ASSAY

SPATIAL DYNAMICS OF
THROMBIN GENERATION
& FIBRIN CLOT GROWTH



SPATIAL DYNAMICS OF COAGULATION

Thrombodynamics is the only laboratory test with adequate physiological model based on the up-to-date understanding of the spatial aspects of in vivo coagulation process.

Thrombodynamics imitates in vitro physiological and pathophysiological processes that occur *in vivo* during hemostatic plug formation or thrombosis.

Unlike other routine coagulation assays the fibrin clot growth process in thrombodynamics assay develops in space and time rather than only in time.

The fibrin clot starts to form, growing from the tissue factor bearing surface, but then propagates into the bulk of the plasma sample without interaction with activator.

Result: 30 min

Distance from the activator

TF-coated Surface

Plasma sample

Tlag [min], Lagtime

V [um/min], Growth rate


CS [um], Clot size

Tsp [min], Spontaneous clots formation time



Thromb Res. 2015: HEMOSTASIS AND THROMBOSIS BEYOND BIOCHEMISTRY: ROLES OF GEOMETRY, FLOW AND DIFFUSION

EFFECT OF ANTICOAGULANTS

Thrombodynamics is sensitive to all types of anticoagulants. Thrombodynamics is more sensitive to heparin than APTT, comparable to anti-Xa and higher than that of TGT and TEG.




J Perinat Med 2017: THE LABORATORY CONTROL OF ANTICOAGULANT THROMBOPROPHYLAXIS DURING THE EARLY POSTPARTUM PERIOD AFTER CESAREAN DELIVERY



COAGULATION STATES AND FIBRINOLYSIS

Thrombodynamics is both a qualitative and quantitative evaluation of the coagulation status and fibrinolysis. Clot formation/lisys process is registered by a CCD camera and then calculated.



LOT [min], Lysis onset time

26,8

LP [%/min], Lysis progretion

3,6

Talanta 2018: THROMBODYNAMICS, A NEW GLOBAL COAGULATION TEST: MEASUREMENT OF HEPARIN EFFICIENCY

WWW.THROMBODYNAMICS.COM

MAIN PARAMETERS AND INTERPRETATION

Tlag, [min] - Lag-time - time between contact of plasma sample with activator and start of clot growth. This parameter is sensitive to the initial stage of blood coagulation and reactions of the extrinsic pathway. Prolongation of this parameter is caused by hypocoagulation of differing nature: deficiency of factors VII and X, (direct thrombin or factor Xa inhibitors, vitamin K antagonists). Shortening is rarely observed, and can be due to different causes of hypercoagulation.

Vi, [um/min] - Initial rate of clot growth - calculated on the interval 2-6 minutes after the beginning of clot growth, also describes initial stages of clot growth but it is spatial elongation rather than local increase of thrombin concentration. A low Vi indicates differing hypocoagulation states (factors VII or X deficiency, anticoagulant agents – factor Xa inhibitors, thrombin inhibitors, vitamin K-antagonists, UFH and LMWH). A high Vi indicates differing hypercoagulation states.

Vst, [um/min] - Rate of clot growth - is the average rate of clot growth calculated on the interval 15-25 min after the beginning of clot growth. If there are no spontaneous clots Vst and V are equal. In the presence of active spontaneous clotting Vst is not calculated. This parameter characterizes the propagation stage of blood coagulation and it is sensitive to all coagulation cascade reactions, including the contact pathway and excluding the initiation reactions of the extrinsic pathway. Decreased Vst indicates various different hypocoagulation states (factors V, VIII, IX, X, XI or thrombin deficiency; anticoagulant agents – vitamin K antagonists, UFH and LMWH). Increased Vst value indicates different hypercoagulation states.

Tsp, [min] - Spontaneous clots - formation time - is the time that spontaneous clots appear in the sample volume which had no initial contact with the activating insert, characterizes clotting independent of the activator surface. Under normal condition, no spontaneous clotting is observed. Spontaneous clotting is induced by circulating activators, active coagulation factors and microparticles. Indicates a high pro-thrombotic tendency.

CS, [um] - Clot size - at the 30th minute of measurement. It is sensitive to all major components and processes of blood coagulation, because it is defined by both Tlag and rate of clot growth.

D, [a.u.] - Clot density - parameter reflects firmness and structure of a formed clot. It reflects quantity and biological activity of fibrinogen, but cannot replace direct measurement of fibrinogen concentration. Sensitive to fibrinogen (concentration and polymerization ability) and factor XIII activity.

Ast, [Activity Unit/L] - Stationary amplitude of thrombin peak - as thrombin generation propagates in space as a moving peak (Dashkevich et al, Biophys J 2012), height of this peak is calculated as a maximal activity of thrombin in the fibrin formation zone which moves from the activator while clot grows. The parameter characterizes the propagation stage of blood coagulation.

Vt, [um/min] - Rate of thrombin peak propagation- characterizes the propagation stage of blood coagulation, sensitive to changes in intrinsic pathway of blood coagulation; factors VIII, IX, XI, V, X and thrombin concentration and also to phospholipid vesicles concentration in plasma.


ETP_ATG, [AU*min/L] - Thrombin potential, **Cmax_ATG, [AU/L]** - Maximum concentration, **Tmax_ATG, [min]**- Time to thrombin peak - are calculated on the activating surface and similar to homogeneous Thrombin Generation Test parameters.

SPATIAL DYNAMICS OF COAGULATION

Thrombodynamics imitates *in vitro* physiological and pathophysiological processes that occur *in vivo* during hemostatic plug formation or thrombosis.

Unlike other routine coagulation assays the fibrin clot growth process in thrombodynamics assay develops in space and time rather than only in time.

The fibrin clot starts to form, growing from the tissue factor bearing surface, but then propagates into the bulk of the plasma sample without interaction with activator.



Result: 30 min

Distance from the activator

TF-coated Surface


Plasma sample

Tlag [min], Lagtime
V [um/min], Growth rate
CS [um], Clot size
Tsp [min], Spontaneous clots formation time

Thromb Res. 2015:

EFFECT OF ANTICOAGULANTS

PLOS ONE 2018: THROMBODYNAMICS—A NEW GLOBAL HEMOSTASIS ASSAY FOR HEPARIN MONITORING IN PATIENTS UNDER THE ANTICOAGULANT TREATMENT



COAGULATION FACTOR DEFICIENCIES

FVIII deficiency (Hemophilia A)

Control

Clot size, μm


Time, min

Spontaneous clots, %

	FVIII deficiency (Hemophilia A)	Control
V	12,9	25,2
Tlag	1,1	0,9
Tsp	None	None
Vi	28,9	48,7
CS	601	1 058
D	17 302	18 159

Biochim Biophys Acta. 2002; 1400: 115-121



HEMOPHILIA A AND B ARE ASSOCIATED WITH ABNORMAL SPATIAL DYNAMICS OF CLOT GROWTH



PHOSPHOLIPID SURFACE AND PLATELETS


The vesicle composition is:
PS:PC:PE = 20:20:60 molar ratio

	4 μM of phospholipids	Control
V	36,5	26,6
Tlag	0,9	0,9
Tsp	None	None
Vi	71,3	58,8
CS	1 525	1 161
D	16 120	20 622



PROCOAGULANT AGENTS

	[TF] = 0,25 pM	Control
V	39,7	34
Tlag	2,4	1,6
Tsp	28	None
Vi	48,8	46,7
CS	1 202	1 119
D	16 616	19 311



Control

Effect of fXla

Clot size, μm

Time, min

Spontaneous clots, %

	Effect of fXla	Control
V	81,4	26
Tlag	0,9	0,9
Tsp	12,4	None
Vi	81,4	51,1
CS	-	1 088
D	27 501	20 622

PLOS ONE 2014;

ACTIVATING MICROPARTICLES TOGETHER WITH FACTORS IXa AND XIa INDUCE SPONTANEOUS CLOTTING IN PLASMA OF HEMATOLOGY AND CARDIOLOGIC PATIENTS

THROMBIN GENERATION

TY AND ROBUSTNESS OF SPATIALLY DEPENDENT THROMBIN GENERATION AND FIBRIN CLOT PROPAGATION

The diagram on the left shows a cross-section of a microfluidic device with two channels. A blue arrow labeled "Fluorescence (thrombin) propagation" points downwards through the channels. To the right, a graph plots "Thrombin, nM" (y-axis, 0 to 600) against "Time, min" (x-axis, 0 to 40). The curve shows a sharp rise to a peak labeled "Cmax_ATG" at approximately 5 minutes, followed by a decay to a baseline labeled "Tmax_ATG" at approximately 10 minutes. The word "SENSITIVITY" is partially visible at the bottom right.

Parameter	CONTROL	UFH 0,04 IU/ml	Dabigatran 70 ng/ml	Rivaroxaban 65 ng/ml
Ast [AU/L]	153,0	31,7	98,4	68,6
Vt [μm/min]	37,0	8,7	28,3	22,9
Tlag [min]	0,8	0,8	1,3	1,8